



# PET life™

## Veterinary Coagulation Analyzer NX-102



### PT/INR, APTT, TT, FIB, ACT

- Sample Type: Citrated Venous Whole Blood or Capillary Blood (PT/INR Only)
- Sample Volume: 20 µL
- Time to Result: ≤5 minutes
- Precision: CV≤5%
- Storage: valid for 18 months at room temperature



### Portable

Small footprint (weighs only 2lbs);  
Uses rechargeable battery

### Convenient

One-step process with no accessories required,  
Maintenance free

### Fast & Reliable

Results in 3 minutes or less  
CV ≤ 5%

### APPLICATIONS

- Screening test to assess the exogenous and endogenous coagulation system function for canine and feline
- Monitoring anticoagulant treatment
- Monitoring thrombolytic therapy
- Monitoring the patients with certain diseases, such as hepatopathy, infection, and burn injury

### Test Items

Cat. No.	Item	Sample Volume	Linear Range	Reference Range Canine	Reference Range Feline	Stored at
VC-PT	PT	20µL	5-90s	5-15s	6-15s	39°~86°F
/	INR	/	/	INR:0.5-1.6	INR:0.6-1.6	/
VC-APTT	APTT	20µL	10-100s	15-45s	15-43s	39°~86°F
VC-FIB	FIB	20µL	1-5g/L	1-3g/L	1-2.5g/L	39°~86°F
VC-ACT	ACT	20µL	50-600s	50-150s	50-160s	39°~86°F
VC-TT	TT	20µL	8-90s	8-20s	9-22s	39°~86°F

\* Whole blood sodium citrate anticoagulation, Note: INR is calculated by PT

※ Analyzers and assays are ordered upon request

Wondfo

Wondfo USA Co., Ltd.

Wondfo USA Co., Ltd.  
Add: 545 Willowbrook Centre Parkway, Unit B, Willowbrook, IL 60527  
Tel: 630-468-2199 Email: info@wondfousa.com Website: www.wondfousa.com