



PETlife™

PETlife Immunofluorescence Analyzer

Rapid Point-of-Care Pet Testing

The PETlife Immunofluorescence Analyzer is designed for convenience and ease-of-use. It is fully calibrated to reduce manual errors and delivers results in minutes. With fully automated touch screen programming, the PETlife Immunofluorescence Analyzer is easy to operate and the small instrument size makes it suitable for use in a variety of settings. Get lab-quality results right where you work.

Easy and Flexible

Simple workflows requiring only ~1 minute hands-on time. Simply add the sample to the test cassette, then insert the cassette into the analyzer and read the result.

Fast and Convenient

No more waiting for lab results; most assay results are available within 15 minutes.

Streamlined Sample Processing

Quick result setting to test multiple pets and traceability to track key data points.



Instrument Overview

- Touchscreen monitor
- ID chip inlet
- Test cassette holder
- Built-in printer
- "Standard" and "Quick Test" modes available

Technical Specifications

- Weight: 6.2lbs/2.8kg
- Dimension: 10.6 x 9.4 x 5.7 inch (270 x 238 x 146mm)
- Voltage input: 100-240V, 50/60Hz, 1.4-0.7A

Software

- Easy to use software program for guided instrument operation, reporting, and troubleshooting
- Self-calibrating system

Rev. A1
Rel. 07/15/2022