



**SAFE life™**

# T-Swab Cotinine (COT) & Alcohol Oral Fluid Test

## What does SAFElife™ T-Swab Cotinine (COT) & Alcohol Oral Fluid Test do?

SAFElife™ T-Swab Cotinine (COT) & Alcohol Oral Fluid Test is a rapid oral fluid test for screening exposure to tobacco smoke including cigarettes and vapor as well as alcohol using.

### COT & Alcohol Detection

Cotinine (COT), a biomarker for tobacco exposure detection, is an alkaloid found in tobacco and is also a major metabolite of Nicotine, which produces stimulation of the autonomic ganglia and central nervous system in humans.

Saliva alcohol test reflect the Blood Alcohol Concentration (BAC), the United States Department of Transportation (DOT) has established a BAC of 0.02% (0.02g/dL) as the cut off level at which an individual is considered positive for the presence of alcohol.

### Optional Alcohol test

Ordering our product with or without Alcohol test upon your demand.  
Optional Alcohol test with transparent window design for clear viewing.

## Choose SAFElife™ T-Swab Cotinine (COT) & Alcohol Oral Fluid Test for the detection of recent Nicotine use and Alcohol use.


- Light weighted pen-shaped device
- No privacy room needed for sample collection
- Simple, rapid, minimally-invasive saliva collection




### Workflow

Ready to use & Soon to get results

Remove the cap




Place the collection pad under the tongue



Get reliable results at 5 minutes for COT test, at 2 minutes for Alcohol test



5 minutes



\* Please refer to the Package Insert for detailed procedure.

Item code	Product name	Description	Complexity	Product image
SODOA-116-COT-P	SAFElife™ T-Swab Cotinine (COT) Oral Fluid Test	1-Panel T-Swab; COT 30 ng/ml	Prescription Use Only	
SODOA-126NT-A	SAFElife™ T-Swab Cotinine (COT) & Alcohol Oral Fluid Test	2-Panel T-Swab; COT+ALC	Forensic Use Only	

Rel: 2024/05/17