



**BABYlife™**

## Pregnancy and Ovulation Tests

BABYlife products provide a simple and easy way to interpret pregnancy and fertility test results with a high degree of accuracy. Our pregnancy and ovulation tests offer different levels of analytical sensitivity. BABYlife products offer standard and early detection, or "higher sensitivity" tests.



### BABYlife Pregnancy Test

*The BABYlife Pregnancy Test is a qualitative test which detects the presence of human chorionic gonadotropin ("HCG"), that appears in urine very early during pregnancy.*

- Results (10mIU HCG: at 5 minutes; 25mIU HCG: in 3-5 minutes; HCG Urine/Serum Combo: in 3-5 minutes.)
- Over 99% Accuracy (from the day of expected period)
- Detects 25 mIU/mL hCG in serum and urine
- Detects 10 mIU/mL hCG in urine
- CLIA waived for urine, Moderate complexity for serum.
- A variety of forms are available: Strip, Cassette, Midstream

### BABYlife Ovulation Test

*The BABYlife Ovulation Test is a qualitative test used to predict when there is a lutenizing hormone (LH) surge, and in turn, when you are likely to ovulate.*

- Results at 5 minutes
- Over 99% accuracy
- Detects 25 mIU/mL LH in urine
- A variety of forms are available: Strip, Cassette, Midstream

Item Code	Description	Box Size / Case Size	Classification
W1-C	hCG Pregnancy Test - Cassette	25 / 20	CLIA WAIVED
W1-S	hCG Pregnancy Test - Strip	20 / 250	CLIA WAIVED
W1-M	hCG Pregnancy Test - Midstream	25 / 20	CLIA WAIVED
W1-M4	hcG Pregnancy Test - Early Result Midstream	25 / 20	CLIA WAIVED
W2-S	LH Ovulation Test - Strip	50 / 100	CLIA WAIVED

Rel.: 2024/08/21