



SAFE life™

Xylazine (XYL) is a non-opioid veterinary tranquilizer, not approved for human use, that is commonly used in veterinary medicine as a sedative and analgesic to calm and restrain animals, particularly horses. However, studies show that people exposed to xylazine often knowingly or unknowingly ingest it in combination with other drugs, particularly illicit fentanyl, which increases the risk of suffering a fatal drug poisoning.

SAFElife™ Xylazine (XYL) Liquid & Powder Drug Test Strip is an immunochromatographic assay for the qualitative determination of Xylazine in liquid, powder or pill samples.



Xylazine (XYL) Liquid & Powder Drug Test Strip


Features and Benefits

- Detects Xylazine in concentrations at 500 ng/mL
- Results in 5 minutes
- Built-in procedural control
- Cost-effective option

Kit Contents

- Individually packaged test strips
- Package Insert

Test Procedure



Item Code	Description	Box Size / Case Size	Complexity
WSPXYL-101	1-Strip; Drug Liquid / Powder Detection; Xylazine 500 ng/ml	25/120	Forensic Use Only