



**SAFE life™**

# C-Cup Multi-Drug Urine Test Cup

## Rapid Point-of-Care Toxicology Testing

- Cost Effective
- Reduce storage space
- Lower shipping expense
- CLIA Waived, 510(k) cleared & Forensic Use Only
- Same technical features as the T-Cup
- Test up to 20 drugs per cup



## CLIA-waived Drugs & Cutoff Levels

|                |                                                          |                 |              |                                      |                 |
|----------------|----------------------------------------------------------|-----------------|--------------|--------------------------------------|-----------------|
| AMP AMP AMP    | Amphetamine (AMP)                                        | 1000, 500 ng/mL | TCA TCA TCA  | Tricyclic Antidepressants (TCA)      | 1000 ng/mL      |
| BAR BAR BAR    | Barbiturates (BAR)                                       | 300 ng/mL       | MET MET MET  | Methamphetamine (MET/mAMP)           | 1000, 500 ng/mL |
| BZO BZO BZO    | Benzodiazepines (BZO)                                    | 300 ng/mL       | MDMA MDMA MD | Methylenedioxymethamphetamine (MDMA) | 500 ng/mL       |
| BUP BUP BUP    | Buprenorphine (BUP)                                      | 10 ng/mL        | MOP MOP MOP  | Morphine (MOP/OPI300)                | 300 ng/mL       |
| COC COC COC    | Cocaine (COC)                                            | 300, 150 ng/mL  | OPI OPI OPI  | Opiate (OPI)                         | 2000 ng/mL      |
| THC THC THC    | Cannabinoids (THC)                                       | 50 ng/mL        | OXY OXY OXY  | Oxycodone (OXY)                      | 100 ng/mL       |
| EDDP EDDP EDDP | 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP) | 300 ng/mL       | PCP PCP PCP  | Phencyclidine (PCP)                  | 25 ng/mL        |
| MTD MTD MTD    | Methadone (MTD)                                          | 300 ng/mL       | PPX PPX PPX  | Propoxyphene (PPX)                   | 300 ng/mL       |

## Forensic Use Panels

|                   |                            |                |             |                                  |                      |
|-------------------|----------------------------|----------------|-------------|----------------------------------|----------------------|
| 6-MAM 6-MAM 6-MAM | Monoacetylmorphine (6-MAM) | 10 ng/mL       | K2 K2 K2    | Synthetic Cannabinoid (K2)       | 50 ng/mL             |
| AMP AMP AMP       | Amphetamine (AMP)          | 300, 500 ng/mL | KET KET KET | Ketamine (KET)                   | 300, 1000 ng/mL      |
| BUP BUP BUP       | Buprenorphine (BUP)        | 5 ng/mL        | KRA KRA KRA | Kratom (KRA)                     | 100, 300 ng/mL       |
| BZO BZO BZO       | Benzodiazepines (BZO)      | 100, 200 ng/mL | LSD LSD LSD | Lysergic Acid Diethylamide (LSD) | 20 ng/mL             |
| COC COC COC       | Cocaine (COC)              | 100, 150 ng/mL | MET MET MET | Methamphetamine (MET/mAMP)       | 300, 500 ng/mL       |
| COT COT COT       | Cotinine (COT)             | 200 ng/mL      | MOP MOP MOP | Morphine (MOP or OPI)            | 100 ng/mL            |
| EDDP EDDP EDDP    | EDDP                       | 100 ng/mL      | MTD MTD MTD | Methadone (MTD)                  | 200 ng/mL            |
| EtG EtG EtG       | Ethyl Glucuronide (EtG)    | 500 ng/mL      | MQL MQL MQL | Methaqualone (MQL)               | 300 ng/mL            |
| FTY FTY FTY       | Fentanyl (FTY)             | 20 ng/mL       | PGB PGB PGB | Pregabalin (PGB)                 | 500 ng/mL            |
| GAB GAB GAB       | Gabapentin (GAB)           | 2000 ng/mL     | THC THC THC | Cannabinoids (THC)               | 15, 25, 40 ng/mL     |
| HMO HMO HMO       | Hydromorphone (HMO)        | 300 ng/mL      | TRA TRA TRA | Tramadol (TRA)                   | 100, 200, 1000 ng/mL |

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



# SAFE life™

## Product Information

※ Kit Size: 25 / Case Size: 4

| Item Code       | Description                                                                                                                | Complexity        |
|-----------------|----------------------------------------------------------------------------------------------------------------------------|-------------------|
| CDOA-254        | 5-Panel C-Cup; AMP, COC, mAMP, OPI, THC                                                                                    | CLIA WAIVED       |
| CDOA-564        | 6-Panel C-Cup; BZO, COC, mAMP, OPI, OXY, THC                                                                               | CLIA WAIVED       |
| CDOA-8105       | 10-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC                                                         | CLIA WAIVED       |
| CDOA-3124       | 12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI, OXY, PCP, TCA, THC                                                | CLIA WAIVED       |
| CDOA-6125       | 12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC                                               | CLIA WAIVED       |
| CDOA-6125A3     | 12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3AD                                      | CLIA WAIVED       |
| CDOA-L6125      | 12-Panel C-Cup, low cutoff; AMP500, BAR, BUP, BZO, COC150, mAMP500, MDMA, MOP, MTD, OXY, PCP, THC                          | CLIA WAIVED       |
| CDOA-9125A3E3F  | 12-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG300, FTY with 3AD                                   | Forensic Use Only |
| CDOA-L9135NTA3  | 13-Panel C-Cup without THC, low cutoff; AMP500, BAR, BUP, BZO, COC150, MDMA, MET500, MOP, MTD, OXY, PCP, PPX, TCA with 3AD | CLIA WAIVED       |
| CDOA-9165EFTK   | 16-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY, TRA, K2                            | Forensic Use Only |
| CDOA-9165EFTKA3 | 16-Panel C-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, OXY, PCP, THC, ETG, FTY, TRA, K2 with 3AD                    | Forensic Use Only |
| CDOA-S1185EFTKK | 18-Panel C-Cup, Special TRA; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC, ETG, FTY, TRA100, K2, KRA  | Forensic Use Only |

Rel: 2024/05/17