


Wondfo®

bringing science to life

2025

USA Product Catalog



Bringing Science to Life. Across America.

Wondfo USA Co., Ltd., Headquartered in Bolingbrook, IL, and founded in 2009, is a leading manufacturer of rapid, diagnostic, over-the-counter and point-of-care tests, and analyzers, for identification of substances and disease conditions, as well as health monitoring. Wondfo USA established a distribution center and office suite in Illinois and an R&D and manufacturing center in California, backed by decades of research experience.

Wondfo USA is a subsidiary of Wondfo Biotech Co., Ltd., which was founded in 1992. In 2015 Wondfo Biotech became the first Chinese point-of-care test manufacturer to be publicly traded, demonstrating the company's wealth of expertise and product innovation.

Wondfo manufactures products in four major categories including, Toxicology Testing (SAFElife™), Infectious Disease (WELLlife™), Women's Reproductive Health (HORMONElife™) and Veterinary Diagnostics (PETlife™). Our products have been distributed globally, benefiting millions of people worldwide.

Wondfo is devoted to the development of innovative products and consistently expanding our portfolio to meet clinical, diagnostic and public health needs. Our mission is to provide solutions that save lives as well as improve day-to-day quality of life.



Wondfo[®]
San Diego

Table of Contents

OEM service is available upon request



WELLlife™

Pages 6 to 11

COVID-19 / Flu A&B Test Solution.....	7
COVID-19 / Flu A&B Home Test.....	8
COVID-19 / Influenza A&B Antigen Test.....	8
COVID-19 Antigen Home Test.....	9
Diagnostic Tests.....	10
Strep A Test.....	11
iFob Test.....	11



SAFElife™

Pages 12 to 31

Toxicology Solutions.....	13
T-Cup Multi-Drug Urine Test Cup.....	15
C-Cup Multi-Drug Urine Test Cup.....	16
F-Cup Multi-Drug Urine Test Cup	17
S-Cup Multi-Drug Urine Test Cup	17
T-Dip Marijuana (THC) Multi-Level Urine Drug Test Panel.....	17
Fentanyl Urine (Home) Test (Cassette/Dip Card)	18
T-Dip Multi-Drug Urine Test Panel.....	19
T-Dip Single Drug Urine Test Card	20
Single Drug Urine Test Strip.....	21
Alcohol Saliva Test	21
T-Square Multi-Drug Oral Fluid Test	22
T-Swab Multi-Drug Oral Fluid Test	23
T-Swab Cotinine (COT) & Alcohol Oral Fluid Test	24
Hair Processor and Multi-Panel Rapid Hair Test.....	25
Harm Reduction Tests.....	26
Xylazine (XYL) Liquid & Powder Drug Test Strip.....	27
Fentanyl (FTY) Liquid & Powder Drug Test Strip	28
Benzodiazepines (BZO) Liquid & Powder Drug Test	29
H1 Portable THC Monitor.....	30
Drug Test Cutoff List	31



HORMONElife™

Pages 32 to 37

Rapid Hormone Tests	33
Pregnancy and Ovulation Tests	34
Ultra Accurate FSH Test Strip	35
PdG Urine Test.....	36
HORMONElife App.....	37



WELLlife™

COVID-19 / Flu A&B Test Solution

This product differentiates between SARS-COV-2, influenza A and influenza B antigens with a single test.

Features & Benefits

- 510(k) cleared for over-the-counter (OTC) use AND Professional use
- Outstanding clinical performance
- Results available in 10 minutes
- Extended detection window
- Wide storage temperature



WELLlife™ COVID-19/Flu A&B Home Test



WELLlife™ COVID-19/Influenza A&B Antigen Test

COVID-19 / Flu A&B Home Test

- Rapid & Easy to Use**

Differentiate between SARS-COV-2, influenza A, and influenza B antigens with a simple nasal swab in just 10 minutes.

- Up-to-Date Variant Detection**

Detects multiple COVID-19 variants, including Omicron and the XEC variant, as well as avian influenza (bird flu).

- Suitable For Ages 2 and Above**

Adult collection is required for testing children 2-13 years old.

- Pre-Loaded 12-Pack Retail Holder Design**

Retail holders come pre-loaded with a total of 12 OTC 1/2ct packs for a total of 12/24 tests. Convenient for retail and pharmacy countertop displays or storage. Features a header card display with QR code for further product details and instructions.

- Extended Detection Window**

Offers a broader detection window within four days of symptom onset.



Item Code	Description	Classification
WV01P0001-1	WELLlife™ COVID-19 / Influenza A&B Home Test 1 Test/Kit	510(k) Cleared for OTC/Home Use
WV01P0001-2	WELLlife™ COVID-19 / Influenza A&B Home Test 2 Tests/Kit	510(k) Cleared for OTC/Home Use
WV01P0001-5	WELLlife™ COVID-19 / Influenza A&B Home Test 5 Tests/Kit	510(k) Cleared for OTC/Home Use
WV01P0001-10	WELLlife™ COVID-19 / Influenza A&B Home Test 10 Tests/Kit	510(k) Cleared for OTC/Home Use
WV01P0001-25	WELLlife™ COVID-19 / Influenza A&B Home Test 25 Tests/Kit	510(k) Cleared for OTC/Home Use

COVID-19 / Influenza A&B Antigen Test

- Rapid Test Results**

Results available in 10 minutes allow for testing and treatment decision-making during the same office visit.

- Up-to-Date Variant Detection**

Detects multiple COVID-19 variants, including Omicron and the XEC variant, as well as avian influenza (bird flu).

- Wide Storage Temperature**

36°F-86°F (2°C-30°C) allows for easier product storage.

- Quality Assurance**

Optional external control testing kit empowers providers to make faster, more confident decisions.

- Extended Detection Window**

Offers a broader detection window within 4 days of symptom onset.



Item Code	Description	Classification
WCOVFLU-POC-5	WELLlife™ COVID-19 / Influenza A&B Antigen Test	510(k) Cleared for Professional Use
WCOVFLU-POC-10	WELLlife™ COVID-19 / Influenza A&B Antigen Test	510(k) Cleared for Professional Use
WCOVFLU-POC-25	WELLlife™ COVID-19 / Influenza A&B Antigen Test	510(k) Cleared for Professional Use
WCOVFLU-CON-5	WELLlife™ COVID-19 / Influenza A&B Antigen Test Control Kit	N/A

COVID-19 Antigen Home Test

510(k) Cleared for Over The Counter (OTC)/Home use

An accurate at-home antigen test provides a fast and reliable way to detect the COVID-19 virus using nasal swabs from individuals showing symptoms. Ideal for quick checks!

Why WELLlife™ COVID-19 Antigen Home Test is the superior choice:

- **Easy to use and rapid results**

A single nasal swab offers accurate results in just 10 minutes, reducing healthcare visits and limiting virus spread, ideal for the elderly and parents with young children.

- **Reliable and accurate**

510(k) cleared confidence, providing the highest levels of performance and accuracy.

- **Pre-Loaded 12-Pack Retail Holder Design**

Retail holders come pre-loaded with a total of 12 OTC 1/2ct packs for a total of 12/24 tests. Convenient for retail and pharmacy countertop displays or storage.



Ordering Information

Item Code	Description	Classification
WL03P0001-1	WELLlife™ COVID-19 Antigen Home Test 1 Test/Kit	510(k) cleared for OTC / Home use
WL03P0001-2	WELLlife™ COVID-19 Antigen Home Test 2 Tests/Kit	510(k) cleared for OTC / Home use
WL03P0001-5	WELLlife™ COVID-19 Antigen Home Test 5 Tests/Kit	510(k) cleared for OTC / Home use
WL03P0001-10	WELLlife™ COVID-19 Antigen Home Test 10 Tests/Kit	510(k) cleared for OTC / Home use
WL03P0001-25	WELLlife™ COVID-19 Antigen Home Test 25 Tests/Kit	510(k) cleared for OTC / Home use

Diagnostic Tests

WELLlife™ provides rapid, point-of-care clinical diagnostic testing products that address critical needs in the market. The CLIA waived WELLlife™ POC diagnostic tests deliver rapid and reliable results, allowing for faster diagnosis and treatment.

Features & Benefits

- Easy to use
- Results available in 10 minutes
- Room temperature storage
- No dietary or medication restrictions
- Easy sample collection



Group A Streptococcus (GAS) Test

Allows providers to rapidly identify the presence of GAS from throat swab specimens and initiate treatment while the patient is still in the office.



Fecal Immunoassay Tests

Aids in the diagnosis of lower gastrointestinal (G.I.) bleeding, as a method of identifying colorectal cancer precursors and other lower G.I. disorders.

Strep A Test

One Step Strep A Swab Test is a rapid, one step lateral flow immunoassay for the qualitative detection of the Group A streptococcus (GAS) bacteria from throat swab specimens.



Item Code	Description	Box Size / Case Size	Classification
PRE-STA-25	Step A Throat Swab Test	25 / 20	CLIA waived

iFOB Immunochemical Fecal Occult Blood Test (FIT)

- *Detects human hemoglobin in stool specimens*
- *No dietary or medication restrictions*
- *Easy sample collection*



Item Code	Description	Box Size / Case Size	Classification
PRE-FIT-25KIT	Kit containing iFOB test cassettes, collection tubes and mailers.	25 / 18	CLIA waived
PRE-FIT-25A	Kit containing iFOB test cassettes and collection tubes.	25 / 20	CLIA waived



SAFElife™

Toxicology Solutions

SAFElife™ toxicology products provide a wide range of tests and configurations to monitor drug use and/or misuse for a variety of settings. Oral fluid tests are used to detect substances within hours of use, while urine tests are used to detect substances within days of use.


Features & Benefits

- Results available in 5 minutes
- Small sample volume requirements
- 510(k) cleared/CLIA waived and Forensic Use Only configurations
- Room temperature storage
- Specific and accurate



Urine Drug Test Cups

Tamper-resistant, available in a range of sizes, shapes and leakproof lids. Customizable panels and optional Vaculids for clean transfer of urine samples. Test for up to 26 drugs. Optional adulteration tests.



T-Dip Urine Drug Tests

Ranging from single substance tests up to 16 different panels. 510(k) cleared / CLIA waived and Forensic Use Only tests are available. Use when you are not concerned about specimen tampering.



Oral Fluid (Saliva) Tests

Full observation during collection allows oral fluid testing to be performed anytime and anywhere, greatly reducing the likelihood of specimen tampering. Test for up to 12 drugs and alcohol.

Toxicology Solutions



Rapid Hair Test Solution

Our on-site, rapid hair testing solution includes 17 drug tests, and boasts an impressive accuracy rate above 95% when compared with Mass Spectrometry testing. No more waiting for negative results – save time and costs by confirming only positive samples. Identify individuals who have been able to "pass" traditional urine tests while misusing substances to get them the help they need.



H1 Portable THC Monitor

Our compact handheld fluorescent device allows for semi-quantitative detection of cannabis concentration in urine. With its user-friendly interface and rapid results, you can conveniently monitor cannabis levels. Take charge of your testing needs with confidence and accuracy using our portable solution.



Harm Reduction Tests

Harm reduction encompasses a range of pragmatic approaches and concepts aimed at mitigating adverse outcomes linked to drug usage. Harm reduction tests promote safety through testing drug samples while upholding the rights and dignity of individuals who engage in drug use.

T-Cup Multi-Drug Urine Test Cup

- Up to 26 drug tests available, customizable
- Results in 5 minutes
- Easy to use
- 24 months shelf life from date of manufacture
- Female friendly opening
- 6 adulteration tests available: CR, GL, NI, OX, pH, SG
(GL, NI and OX available by special request)



※ Kit Size: 25 / Case Size: 4

Item Code	Description	Classification
TDOA-154	5-Panel T-Cup; AMP, COC, OPI, PCP, THC	CLIA waived
TDOA-254	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC	CLIA waived
TDOA-254A3	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC with 3 AD	CLIA waived
TDOA-164	6-Panel T-Cup; AMP, COC, MET, OPI, PCP, THC	CLIA waived
TDOA-264	6-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, THC	CLIA waived
TDOA-274	7-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, OXY, THC	CLIA waived
TDOA-3104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC	CLIA waived
TDOA-3104A3	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC with 3 AD	CLIA waived
TDOA-4104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA waived
TDOA-7104	10-Panel T-Cup; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC	CLIA waived
TDOA-8104	10-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, THC	CLIA waived
TDOA-4124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA waived
TDOA-6124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA waived
TDOA-6124A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC with 3 AD	CLIA waived
TDOA-6125	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA waived
TDOA-6125A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3 AD	CLIA waived
TDOA-6125V	12-Panel T-Cup w/Vacu Lid; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA waived
TDOA-7125	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA waived
TDOA-2135	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA waived
TDOA-2135A3	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA waived
TDOA-6135F	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, FTY	Forensic Use Only
TDOA-7135EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA	Forensic Use Only
TDOA-7135A3EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, MOP, OXY, THC, ETG, FTY, TRA with 3 AD	Forensic Use Only
TDOA-8135E	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG	Forensic Use Only
TDOA-1144	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	CLIA waived
TDOA-1144A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA waived
TDOA-1145A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA waived
TDOA-8145EF	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only
TDOA-9145EFTK	14-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA, K2	Forensic Use Only
TDOAC-L1223A5EFTKK6	22-Panel T-Cup Clear; AMP300, BAR300, BUP10, BZO200, COC100, EDDP300, HMO300, KET1000, MDMA500, MET300, MOP100, MTD200, OXY100, PCP25, PPX300, THC25, ETG500, FTY20, TRA200, K2 50, KRA300, 6-MAM10 + (CR, NI, OX, pH, S.G.)	Forensic Use Only

C-Cup Multi-Drug Urine Test Cup

- Up to 20 drug tests available, customizable
- Cost effective
- Reduced storage space requirements
- Lower shipping costs
- 510(k) cleared/CLIA waived and Forensic Use Only
- Same technical features as T-Cup



※ Kit Size: 25 / Case Size: 4

Item Code	Description	Classification
CDOA-254	5-Panel C-Cup; AMP, COC, mAMP, OPI, THC	CLIA waived
CDOA-564	6-Panel C-Cup; BZO, COC, mAMP, OPI, OXY, THC	CLIA waived
CDOA-8105	10-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC	CLIA waived
CDOA-8104	10-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, THC	CLIA waived
CDOA-3104A3	10-Panel C-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC with CR/pH/S.G.	CLIA waived
CDOA-3124	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI, OXY, PCP, TCA, THC	CLIA waived
CDOA-6125	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA waived
CDOA-6125A3	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3AD	CLIA waived
CDOA-L6125	12-Panel C-Cup, low cutoff; AMP500, BAR, BUP, BZO, COC150, mAMP500, MDMA, MOP, MTD, OXY, PCP, THC	CLIA waived
CDOA-9125A3E3F	12-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG300, FTY with 3AD	Forensic Use Only
CDOA-6124	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA waived
CDOA-6124A3	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC with CR/pH/S.G.	CLIA waived
CDOA-L9135NTA3	13-Panel C-Cup without THC, low cutoff; AMP500, BAR, BUP, BZO, COC150, MDMA, MET500, MOP, MTD, OXY, PCP, PPX, TCA with 3AD	CLIA waived
CDOA-2135	13-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA waived
CDOA-2135A3	13-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with CR/pH/S.G.	CLIA waived
CDOA-6135F	13-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, FTY	Forensic Use Only
CDOA-7135EFT	13-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, MOP, OXY, THC, EtG, FTY, TRA	Forensic Use Only
CDOA-1145A3	14-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC with CR/pH/S.G.	CLIA waived
CDOA-9145EFTK	14-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, MOP, OXY, THC, EtG, FTY, TRA, K2	Forensic Use Only
CDOA-8145EF	14-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, EtG, FTY	Forensic Use Only
CDOA-1144	14-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	CLIA waived
CDOA-1144A3	14-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC with CR/pH/S.G.	CLIA waived
CDOA-9165EFTK	16-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY, TRA, K2	Forensic Use Only
CDOA-9165EFTKA3	16-Panel C-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, OXY, PCP, THC, ETG, FTY, TRA, K2 with 3AD	Forensic Use Only
CDOA-S1185EFTKK	18-Panel C-Cup, Special TRA; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC, ETG, FTY, TRA100, K2, KRA	Forensic Use Only

F-Cup Multi-Drug Urine Test Cup

- Up to 16 drug tests available, customizable.
- Two flat sides allowing for photocopying of results
- 510(k) cleared / CLIA waived
- Same technical features as T-Cup

※ Kit Size: 25 / Case Size: 4



Item Code	Description	Classification
FDOA-6125	12-Panel F-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC	CLIA waived

S-Cup Multi-Drug Urine Test Cup

- Four flat sides allowing for photocopying of results
- Up to 18 drug tests available, customizable
- Customizable panels
- 510(k) cleared / CLIA waived and Forensic Use Only options
- Same technical features as T-Cup



T-Dip Marijuana (THC) Multi-Level Urine Drug Test Panel

- Semi-quantitative THC detection
- Narrow down the level of THC in a urine sample
- 5 cutoffs available
- Result in 5 minutes

※ Kit Size: 25 / Case Size: 20



Item Code	Description	Classification
WDTH-134M	3-Panel Dip Card; 15/50/100 ng/mL THC cutoffs	Forensic Use Only
WDTH-151	5-Panel Dip Card; 15/50/100/200/300 ng/mL THC cutoffs	Forensic Use Only

Fentanyl Urine Test (Cassette/Dip Card)

For OTC (Home) use and Professional use.

The necessity of ultra-sensitive fentanyl detection

Fentanyl is 50 times more potent than heroin and 100 times more potent than morphine. 2mg of fentanyl, which is 10 times less than the weight of a single grain of rice, is considered a lethal dose to an average adult.¹

For OTC Use



Fentanyl Urine Home Test (Cassette)



Fentanyl Urine Home Test (Dip Card)

For Professional Use, CLIA waived



Fentanyl Urine Test (Cassette), 25 Tests /Kit



Fentanyl Urine Test (Dip Card), 25 Tests /Kit

Product Highlights

- Ultra-Sensitive Fentanyl Urine Detection at 1 ng/mL
- Simple to use, results are a few steps away
- Over 95% accuracy when compared with the gold standard LCMS, ensuring reliable results.
- 510(k) cleared for confidence


Item Code	Description	Kit Size	Classification
WFTY-OTC1-1BX	Fentanyl 1ng Dip Card	1 Test	OTC
WDFT-201-POC1W-25BX	Fentanyl 1ng Dip Card	25 Tests	CLIA waived
WDFT-201-OTC1-1BX	Fentanyl 1ng Cassette	1 Test	OTC
WFTY-POC1W-25BX	Fentanyl 1ng Cassette	25 Tests	CLIA waived

Rel.: 2025/06/26

[1] Deadly Dose, <https://www.dea.gov/onepill>

T-Dip Multi-Drug Urine Test Panel

- Up to 16 drug tests available, customizable
- Various formats and combinations
- Results in 5 minutes
- Room temperature storage
- 24 months shelf life
- Optional built-in adulteration indicator for OEM products



✂ Kit Size: 25 / Case Size: 20

Item Code	Description	Classification
WDOA-234	3-Panel Dip Card; COC,mAMP,THC	CLIA waived
WDOA-154	5-Panel Dip Card; AMP,COC,OPI,PCP,THC	CLIA waived
WDOA-254	5-Panel Dip Card; AMP,COC,mAMP,OPI,THC	CLIA waived
WDOA-654	5-Panel Dip Card; BZO,COC,mAMP,OPI,THC	CLIA waived
WDOA-754	5-Panel Dip Card; AMP,BZO,COC,OPI,THC	CLIA waived
WDOA-264	6-Panel Dip Card; AMP,BZO,COC,mAMP,OPI,THC	CLIA waived
WDOA-564	6-Panel Dip Card; BZO,COC,mAMP,OPI,OXY,THC	CLIA waived
WDOA-865NT	6-Panel Card without THC; BUP,BZO,COC,MOP,MTD,OXY	CLIA waived
WDOA-274	7-Panel Dip Card; AMP,BZO,COC,mAMP,OPI,OXY,THC	CLIA waived
WDOA-295	9-Panel Dip Card; AMP,BUP,BZO,COC,MOP,MTD,OXY,PPX,THC	CLIA waived
WDOA-994NT	9-Panel Card without THC; AMP,BAR,BUP,BZO,COC,MET,MTD,OPI,OXY	CLIA waived
WDOA-1104	10-Panel Dip Card; AMP,BAR,BZO,COC,mAMP,MTD,OPI,PCP,TCA,THC	CLIA waived
WDOA-3104	10-Panel Dip Card; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,PCP,THC	CLIA waived
WDOA-4104	10-Panel Dip Card; AMP,BAR,BZO,COC,mAMP,MTD,OPI,OXY,PCP,THC	CLIA waived
WDOA-7104	10-Panel Dip Card; BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,TCA,THC	CLIA waived
WDOA-9104	10-Panel Dip Card; AMP,BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,THC	CLIA waived
WDOA-1124	12-Panel Dip Card; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,PPX,THC	CLIA waived
WDOA-3124	12-Panel Dip Card; AMP,BAR,BUP,BZO,COC,MDMA,MTD,TCA,OPI,OXY,PCP,THC	CLIA waived
WDOA-6124	12-Panel Dip Card; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,THC	CLIA waived
WDOA-6125	12-Panel Dip Card; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC	CLIA waived
WDOA-7125	12-Panel Dip Card; AMP,BAR,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC	CLIA waived
WDOA-814SEF	14-Panel Dip Card; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only

T-Dip Single Drug Urine Test Card

- **Results in 5 minutes**
- **Customizable cut-off values available**
- **Room temperature storage**
- **24 months shelf life**




※ **Kit Size: 25 / Case Size: 20**

Item Code	Description	Classification
WDBU-114	1-Panel Dip Card; BUP Buprenorphine	CLIA waived
WDBZ-114	1-Panel Dip Card; BZO Benzodiazepines	CLIA waived
WDCO-114	1-Panel Dip Card; COC Cocaine	CLIA waived
WDMA-114	1-Panel Dip Card; mAMP/MET Methamphetamine	CLIA waived
WDMT-114	1-Panel Dip Card; MTD Methadone	CLIA waived
WDOP-114	1-Panel Dip Card; OPI Opiates 2000 ng/mL	CLIA waived
WDOX-114	1-Panel Dip Card; OXY Oxycodone	CLIA waived
WDPC-114	1-Panel Dip Card; PCP Phencyclidine	CLIA waived
WDTH-114	1-Panel Dip Card; THC Marijuana/Cannabinoids	CLIA waived
W6AM-114	1-Panel Dip Card; 6-MAM Monoacetylmorphine	Forensic Use Only
WCOT-114	1-Panel Dip Card; COT Cotinine	Forensic Use Only
WDK2-114	1-Panel Dip Card; K2 Synthetic Cannabinoid	Forensic Use Only
** WDTH-115	1-Panel Dip Card; THC Marijuana 15 ng/mL	Forensic Use Only
WDTR-114	1-Panel Dip Card; TRA Tramadol	Forensic Use Only
WETG-114	1-Panel Dip Card; EtG Ethyl Glucuronide	Forensic Use Only
WETG-300	1-Panel Dip Card; EtG Ethyl Glucuronide 300 ng/mL	Forensic Use Only
WFTY-114	1-Panel Dip Card; FTY Fentanyl	Forensic Use Only
** WGAB-114	1-Panel Dip Card; GAB Gabapentin	Forensic Use Only
** WHMO-114	1-Panel Dip Card; HMO Hydromorphone	Forensic Use Only
WKET-114	1-Panel Dip Card; KET Ketamine 1000ng/mL	Forensic Use Only
WKRA-114	1-Panel Dip Card; KRA Kratom	Forensic Use Only
WMQL-114	1-Panel Dip Card; MQL Methaqualone	Forensic Use Only
** WPGB-114	1-Panel Dip Card; PGB Pregabalin	Forensic Use Only

** Ordered upon request

Single Drug Urine Test Strip



※ Kit Size: 100 / Case Size: 50

Item Code	Description	Classification
WDTH-101	1-Strip; THC Marijuana / Cannabinoids	CLIA waived

Alcohol Saliva Test

One Step Alcohol Saliva Test Strip is intended for use as a rapid method to detect the presence of alcohol in saliva, equating to a blood alcohol concentration (BAC) greater than 0.02%.

To confirm the concentration of a positive specimen, an alternate, non-enzymatic technology such as headspace gas chromatography should be used.



Item Code	Description	Tests Per Box / Case Size	Classification
W53-S-100BX	Alcohol Saliva Test	100 / 25	OTC
W53-S-25BX	Alcohol Saliva Test	25 / 168	OTC

T-Square Multi-Drug Oral Fluid Test

- Available for both Forensic Use Only and Employment and Insurance (E&I) Testing
- Optional saturation indicator to confirm sample adequacy
- Fast and easy collection; No privacy required
- Up to 12 different drugs and testing panels including alcohol
- Useful in detecting substance use during the window of impairment



※ Kit Size: 25 / Case Size: 4

Item Code	Description	Classification
QODOA-156NTEUO-I	5-Panel T-Square no THC with indicator; AMP, COC, MET, OPI, PCP	E&I
QODOA-256EUO-I	5-Panel T-Square with indicator; AMP, COC, MET, OPI, THC	E&I
QODOA-166EUO-I	6-Panel T-Square with indicator; AMP, COC, MET, OPI, PCP, THC	E&I
QODOA-366NTEUO-I	6-Panel T-Square no THC with Indicator; AMP, COC, MET, OPI, OXY, PCP	E&I
QODOA-666EUO-I	6-Panel T-Square with indicator; AMP, COC, MET, OPI, OXY, THC	E&I
QODOA-376EUO-I	7-Panel T-Square with Indicator; AMP, COC, MET, OPI, OXY, PCP, THC	E&I
QODOA-496NTEUO-I	9-Panel T-Square no THC with Indicator; AMP, BAR, COC, MDMA, MET, MTD, OPI, OXY, PCP	E&I
QODOA-10106EUO-I	10-Panel T-Square with indicator; AMP, BAR, COC, MDMA, MET, MTD, OPI, OXY, PCP, THC	E&I
QODOA-266-I	6-Panel T-Square with indicator; AMP, BZO, COC, MET, OPI, THC	Forensic Use Only
QODOA-5106	10-Panel T-Square; AMP, BAR, BUP, BZO, COC, MET, OPI, OXY, PCP, THC	Forensic Use Only
QODOA-3126-A*	12-Panel T-Square; AMP, BAR, BUP, BZO, COC, MET, MTD, OPI, OXY, PCP, THC+ALC	Forensic Use Only
QODOA-6126-I	12-Panel T-Square with Indicator; AMP, BAR, BUP, BZO, COC, MDMA, MET, MTD, OPI, OXY, PCP, THC	Forensic Use Only

* "A" at the end stands for alcohol

※ Customizable drug panels & Alcohol test are available.

T-Swab Multi-Drug Oral Fluid Test

Experience the ultimate convenience with our lightweight and compact device.

- **Lightweight & compact**
- **Easy to use**
- **Tests can be run within 2 minutes**
- **Result in 5 minutes**
- **Screens for 7 common drugs of misuse**
- **For Employment and Insurance (E&I) Testing**

Workplace Drug Testing Situations:

Follow-Up, For Cause with suspicion, periodic, post-accident, pre-employment, random, return to duty. Not applicable for Federally regulated workplace and DOT testing programs.

T-Swab Cutoff Levels

Drug Name	Cutoff (ng/mL)
AMP (Amphetamines)	50
COC (Cocaine)	20
THC (Marijuana)	40
MET (Methamphetamine)	50
OPI (Opiates)	40
OXY (Oxycodone)	20
PCP (Phencyclidine)	10

Item Code	Test Combination
** SODOA-156NTEUO	AMP-COC-OPI; MET-PCP
** SODOA-256EUO	AMP-COC-OPI; MET-THC
SODOA-166EUO	AMP-COC-OPI; MET-THC-PCP
** SODOA-666EUO	AMP-COC-OPI; MET-THC-OXY

** Special Order Items



T-Swab Cotinine & Alcohol Oral Fluid Test

What does SAFElife™ T-Swab Cotinine (COT) & Alcohol (ALC) Oral Fluid Test do?

SAFElife™ T-Swab Cotinine (COT) & Alcohol Oral Fluid Test is a rapid oral fluid screening test for tobacco use including cigarettes and vapor as well as alcohol use.

Cotinine & Alcohol Detection

Cotinine (COT), a biomarker for tobacco use, is an alkaloid found in tobacco and is also a major metabolite of Nicotine, which produces stimulation of the autonomic ganglia and central nervous system in humans.

The Saliva alcohol test reflects a Blood Alcohol Contraction (BAC of 0.02% (0.02 g/dL)). This is equivalent to the level at which an individual is considered positive for the presence of alcohol by the United States Department of Transportation (DOT).

Optional Alcohol test

Choose a product with COT and ALC combined or COT only according to your needs.

Product Highlights

- Light weight, flat device
- No privacy needed for sample collection
- Simple, rapid, non-invasive saliva collection
- Results available in 5 minutes



Item code	Description	Classification
** SODOA-116-COT-P	1-Panel COT T-Swab 30 ng/mL	Prescription Use Only
** SODOA-126NT-A	2-Panel COT + ALC T-Swab 30 ng/mL, ALC 0.02%	Forensic Use Only

** Special Order Items

Hair Processor and Multi-Panel Rapid Hair Test

- *Rapid detection of substance use between 14 and 90 days depending on the length of the hair sample.*
- *E&I and Forensic Use Only testing panels*
- *Choose from a comprehensive menu of 18 accurate and reliable drug tests*
- *Customizable panel options*
- *Private and observed sample collection for donor protection*
- *Minimizes the risk of adulteration, ensuring result reliability*
- *Results in 5 mins (15 minutes for THC test)*



Maximum panel size: 12-panel test.

- | | | |
|--------------------------------|---|------------------------|
| • 6-Monoacetylmorphine (6-MAM) | • Cocaine (COC)* | • Methadone (MTD) |
| • Amphetamine (AMP)* | • Fentanyl (FTY)* | • Opiate (OPI)* |
| • Barbitol (BAR) | • Synthetic Marijuana (K2) | • Oxycodone (OXY)* |
| • Benzoylcegonine (BE)* | • Ketamine (KET) | • Phencyclidine (PCP)* |
| • Buprenorphine (BUP) | • Methamphetamine (MET/mAMP200)* | • Marijuana 200 (THC) |
| • Benzodiazepines (BZO)* | • Methamphetamine (MET/mAMP500)* | • Marijuana 100 (THC) |
| | • Methylenedioxymethamphetamine (MDMA)* | • Tramadol (TRA) |


* E&I, others are FUO

Workflow


Efficient, accurate, on-site simplicity for drug testing confidence.




Collect the hair sample, place it into the processing tube, and screw the cap.



Insert the processing tube into the SAFE life™ hair processor and press the start button.



Remove the processing tube. Wait 3 minutes for sample cooling and pour the tube contents into the testing cube.



Read the results between 5 to 15 minutes depending on the tests chosen.

※ **Kit Size: 25 / Case Size: 4**


Item Code	Description	Classification
QHDOA-118-THC	1-Panel Safelife Hair Test; THC 100	Forensic Use Only
QHDOA-128-THC	2-Panel Safelife Hair Test THC 100 and THC 200 cut-off	Forensic Use Only
QHDOA-198	9-Panel SAFelife Hair Test; AMP500, BZO, COC, FTY, MET500, MDMA, OPI, OXY, PCP	Forensic Use Only
QHDOA-1108	10-Panel Safelife Hair Test; AMP, BAR, BUP, BZO, COC, FTY, MAMP, OPI, OXY, PCP	Forensic Use Only
QHDOA-1128	12-Panel Safelife Hair Test; 6-MAM, AMP, BAR, BZO, COC, FTY, mAMP, MDMA, MTD, OPI, OXY, PCP	Forensic Use Only

Harm Reduction Tests

Harm reduction is an evidence-based approach that is critical to engaging with individuals who use drugs, and equipping them with life-saving tools and information to create positive change in their lives and protect them from harm. Harm reduction is a key pillar in the U.S. Department of Health and Human Services' Overdose Prevention Strategy.


Source: SAMHSA (Substance Abuse and Mental Health Services Administration)

SAFElife™ Harm Reduction Liquid & Powder Tests can determine if drugs are present in liquid or powder samples, empowering users with vital decision making information to assist with harm reduction.




Xylazine (XYL) Liquid & Powder Drug Test

Xylazine (XYL) is a non-opioid veterinary tranquilizer, not approved for human use. It has been increasingly linked to overdose deaths. Drug onset is rapid with effects lasting for hours. Studies have shown that individuals knowingly or unknowingly ingest it in combination with other drugs, including fentanyl, commonly resulting in multi-drug toxicity.



Fentanyl (FTY) Liquid & Powder Drug Test

Fentanyl (FTY) is a potent, synthetic opioid analgesic with a rapid onset and short duration of action. Fentanyl is approximately 50 times stronger than heroin and 100 times stronger than morphine. It has become a major contributor to fatal and non fatal overdoses in the U.S.



Benzodiazepines (BZO) Liquid & Powder Drug Test


Benzodiazepines (BZO) are widely used in the United States for their medical benefits, including treating anxiety, insomnia, seizures, and muscle spasms. However, their potential for misuse and dependence has led to a growing focus on harm reduction strategies.

Xylazine (XYL) Liquid & Powder Drug Test


Choose the Test Kit (test strip and mixing accessories) or Test Strip according to your specific requirements.

- Detects Xylazine at concentrations at 500 ng/mL or higher.
- Results as fast as one minute
- Easy to use: Step by step instructions and clearly marked kit components
- Cost-effective: Different kit options to meet diverse needs.
- Forensic Use Only/Harm Reduction Use

All-In-One Kit



Test Kit Pro



Test Kit Premier



Test Strip


Item Code	Description	Classification
WSPXYL-101	1 Xylazine, 500ng/mL test strip for liquid or powder detection	Forensic Use Only, Harm Reduction Use
WSPXYL-101-1PAP	Pro kit with 1 xylazine, 500ng/mL test strip and mixing accessories for drug surface/powder/pill testing	Forensic Use Only, Harm Reduction Use
WSPXYL-101-25PA	Premier kit with 1 xylazine, 500ng/mL test strip, mixing vial and accessories for drug surface/powder/pill testing	Forensic Use Only, Harm Reduction Use

Fentanyl (FTY) Liquid & Powder Drug Test


Choose the Test Kit (test strip and mixing accessories) or Test Strip according to your specific requirements.

- Detects multiple fentanyl analogues in concentrations at 10 ng/mL or higher
- No false positive with MDMA, Methamphetamine, Cocaine or Morphine
- Easy to use: Step by step instructions and clearly marked kit components
- Results as fast as one minute
- Cost-effective: Different kit options to meet diverse needs
- Forensic Use Only/Harm Reduction Use

All-In-One Kit



Test Kit Pro



Test Kit Premier



Test Strip


Item Code	Description	Classification
WSPTY-101	1 Fentanyl, 10ng/mL test strip for liquid or powder detection	Forensic Use Only, Harm Reduction Use
WSPTY-101-1PAP	Pro kit with 1 Fentanyl, 10ng/mL test strip and mixing accessories for drug surface/powder/pill testing	Forensic Use Only, Harm Reduction Use

Benzodiazepines (BZO) Liquid & Powder Drug Test

Choose the Test Kit (test strip and mixing accessories) or Test Strip according to your specific requirements.

- Detects multiple Benzodiazepine analogs
- Easy to use: Step by step instructions and clearly marked kit components
- Results as fast as one minute
- Cost-effective: Different kit options to meet diverse needs
- Forensic Use Only/Harm Reduction Use

All-In-One Kit



Item Code	Description	Classification
WSPBZO-101	1 Benzodiazepines test strip for liquid or powder detection	Forensic Use Only, Harm Reduction Use
WSBZO-101-1PAP	Pro kit with 1 Benzodiazepines test strip and mixing accessories for drug surface/powder/pill testing	Forensic Use Only, Harm Reduction Use

H1 Portable THC Monitor

- *Determine THC concentration in urine from a range of six levels*
- *Tracks the level of THC in urine from hours to days after consumption.*
- *Ultra-accurate results based on laboratory testing technology.*
- *Over 99% performance agreement with clinical data*
- *Results in 5 minutes*

Instrument Parameters

Name	H1 Portable THC Monitor
Analyzer dimensions	4.5 x 2.75 x 2 inches
Screen size	1.5 Inches
Rechargeable battery	24 weeks of standby power.
Weight	≤0.5 lb (including built-in battery)
Data storage	Internal memory stores results for later download
Operating temperature	59°F~86°F (15°C~30°C)

Used in conjunction with the SAFElife™ THC Urine Test Cassette.

Clear interpretation

Level	THC Concentration Detected
L1	0-15 ng/mL
L2	15-50 ng/mL
L3	50-100 ng/mL
L4	100-200 ng/mL
L5	200-300 ng/mL
L6	>300 ng/mL



Item code	Description	Kit size	Format
W947	SAFElife™ H1 Portable THC Monitor	1	Analyzer
W947-THC161	THC Urine Test Cassette for SAFElife™ H1 Portable THC Monitor;	5 Tests	Cassette

Drug Test Cutoff List (USA)

Drug	Abbreviation	Urine - Regular		Urine - Special			Oral Fluid			Harm reduction
		CLIA waived	Forensic Use	CLIA waived	510 (k) cleared	Forensic Use	EGL (Employment & Insurance)	Forensic Regular	Forensic Special	
Amphetamines	AMP	1000		500	300		50	50	40	
Barbiturates	BAR	300					20	60		
Benzodiazepines	BZO	300			200	100 ♦		30		
Buprenorphine	BUP	10				5		5		
Carisoprodol	CAR/SOMA					2000 NEW				
Cocaine	COC	300		150	100		20	20		
Cotinine	COT		200		200 ●					
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	EDDP	300				300				
Ethyl Glucuronide	EtG		500							
Etomidate	ETO					100 NEW				
Fentanyl	FTY		20 ▲	1 ■				30		10
Gabapentin	GAB		2000							
Hydromorphone	HMO		300							
γ-Hydroxybutyric acid	GHB					10µg/mL NEW				
Ketamine	KET					1000				
Kratom	KRA		300			100				
Lysergic Acid Diethylamide	LSD		20							
Marijuana	THC	50 ★			40 ★	15 & 25 ★	40 ★★	25 ★		
Methadone	MTD	300			200		30	30		
Methamphetamine	MET/mAMP	1000		500	300 ♦		50	50		
Methaqualone	MQL		300							
Methcathinone	MCAT					1000 NEW				
Methylenedioxymethamphetamine	MDMA	500					50	100		
Monoacetylmorphine	6-MAM		10							
Morphine	MOP/OPI300	300			100					
Opiates	OPI	2000					40	40	15	
Oxycodone	OXY	100					20	20		
Phencyclidine	PCP	25					10	10		
Pregabalin	PGB		500							
Propoxyphene	PPX	300								
Psilocybin	PSY					2000 NEW				
Synthetic Cannabinoids	K2		50							
Synthetic Cannabinoids	K3					50 NEW				
Tramadol	TRA		200			100 & 1000 ♦				
Tricyclic Antidepressants	TCA	1000								
Xylazine	XYL		500 NEW							500
Alcohol	ALC							0.02 g/dL		
Alcohol	EtOH					0.04 g/dL				

Adulteration test options: Creatinine (CR), Glutaraldehyde (GL), Specific Gravity (S.G.), pH, Oxidant (OX) & Nitrite (NI)

★ 11-nor-Δ⁹-THC-9-COOH

★★ Δ⁹-THC

● Class I exempt; Prescription Use Only

▲ Norfentanyl

■ Fentanyl

♦ Please contact us prior to ordering





HORMONELife™

Rapid Hormone Tests

HORMONElife™ products provide a simple and easy way to interpret pregnancy and fertility test results with a high degree of accuracy. Our pregnancy and ovulation tests offer different levels of analytical sensitivity with standard and early detection, or “higher sensitivity” options.

Features & Benefits

- Results available in 3-5 minutes
- Over 99% accuracy
(from the first day of a missed period)
- Standard and early detection
- Various forms available: strips, cassette, midstream



Pregnancy Tests




Ovulation Test



FSH Test Strip



PdG Urine Test



HORMONElife™ APP

Pregnancy and Ovulation Tests



HORMONElife™ Pregnancy Test

The HORMONElife™ Pregnancy Test is a qualitative test which detects the presence of human chorionic gonadotropin ("hCG"), hormone which appears in urine very early during pregnancy.

- Detects 10 mIU/mL hCG in urine
- Detects 25 mIU/mL hCG in serum and urine
- Over 99% accuracy (from the first day of a missed period)
- Results available within 5 minutes
- CLIA waived for urine, Moderate Classification for serum
- Various forms available: Strip, Cassette, Midstream

HORMONElife™ Ovulation Test

The HORMONElife™ Ovulation Test is a qualitative test used to predict an increase in lutenizing hormone (LH), and in turn, the likely time of ovulation

- Detects 25 mIU/mL LH in urine
- Over 99% accuracy
- Results available in 5 minutes
- 510(k) exempt
- Various forms available: Strip, Cassette, Midstream

Item Code	Description	Kit Size / Case Size	Classification
W1-C-CLIA	1-Cassette; hCG Pregnancy Test	25 / 20	CLIA waived
W1-S	1-Strip; hCG Pregnancy Test	20 / 250	OTC
W2-S	1-Strip; LH Ovulation Test	50 / 100	OTC
W1-M	1-Midstream; hCG Pregnancy Test	25 / 20	OTC
W1-M4	1-Midstream; hCG Pregnancy Test	25 / 20	OTC

Ultra Accurate FSH Test Strip

The HORMONElife™ Ultra Accurate FSH Test Strip is a rapid chromatographic immunoassay for detection of Follicle Stimulating Hormone (FSH) in urine. The FSH test helps women determine their fertility and menopausal status, key indicators for women attempting to conceive.

- **Rapid and Easy to Use**

Results in just 10 minutes with a simple urine dip test.

- **High Sensitivity**

Semi-quantitative results with varying sensitivity levels to cover different stages of the menstrual cycle: 5/10/25/50 mIU/mL.

- **Health Empowered**

Gain valuable insights into fertility status.

- **510(k) Exempt**

Available without prescription.



Item Code	Description	Kit Size / Case Size	Classification
W52-S-SQ-25BX	Semi-Quantitative Follicle Stimulating Hormone (FSH) Strips 5/10/25/50 mIU/mL	25 Tests / 20 Boxes per Case	OTC

PdG Urine Test

Pregnanediol Glucuronide (PdG) is the major urine metabolite of progesterone. The HORMONElife™ PdG Urine Test Strip detects elevated levels of PdG in urine. Using the test at various specific intervals helps confirm ovulation.

Ovulation & Fertility

Pregnancy cannot occur without ovulation, but ovulation doesn't necessarily occur every cycle for every woman. It's estimated that in a third of all menstrual cycles, ovulation doesn't occur.^[1]


Women struggling to conceive need to understand their fertility and ovulation status.

PdG test can provide insights into progesterone levels and overall reproductive health.

What is PdG?

Pregnanediol Glucuronide (PdG) is the major urine metabolite of progesterone, that helps confirm whether ovulation has occurred. Studies show that PdG levels in urine correlate with progesterone levels in blood.^[2]

Progesterone and PdG levels typically peak around 5 to 9 days after ovulation. This is the peak fertility window. It may take up to 5 days after peak fertility for PdG to return to a sufficient level at which to produce a positive test result.



Item Code	Description	Kit Size / Case Size	Classification
DW199-S-6BX-A	PdG Urine Test	6 Tests / 160 Boxes per Case	OTC
DW199-S-12BX-A	PdG Urine Test	12 Tests / 204 Boxes per Case	OTC

Features & Benefits

- Helps to Confirm Ovulation
- High Sensitivity
- 510(K) Exempt
- Rapid & Easy to Use
- Results in just 5 minutes with a simple dip test

Rel.: 2025/06/26

[1] Prior, Jerilynn C et al. Ovulation Prevalence in Women with Spontaneous Normal-Length Menstrual Cycles - A Population-Based Cohort from HUNT3, Norway, 2015



[2] Monitoring the menstrual cycle: Comparison of urinary and serum reproductive hormones referenced to true ovulation.

HORMONElife™ App


Period, Ovulation, and Perimenopause Tracker

HORMONElife™ is the ultimate fertility tracking app that empowers women on the journey to parenthood, with unmatched precision and simplicity.


- **No More Guessing**
Automatically scan ovulation tests for a clear digital interpretation.
- **Track Hormone Trends**
Monitor hormone levels with intuitive charts to understand patterns.
- **Predict Fertile Window**
Intelligent algorithms forecast the optimal conception window with precision.
- **Boost Conception Chances**
Accurate predictions based on hormone levels, temperature, and age to enhance your likelihood of conception.



Capture LH, FSH & PdG result levels with a phone.



Track hormone changes with charts.



Track period and ovulation days.



Wondfo[®] || bringing science to life

+1 (630) 468-2199 || wondfo@wondfousa.com || **WondfoUSA.com**